

Hirudoid® forte

Cream

Mucopolysaccharide polysulphate

Antithrombotic and anti-inflammatory preparation for cutaneous treatment of

inflammation of the veins such as thrombophlebitis, haematoma and for softening of hard scars.

Hirudoid forte cream contains Organo-heparinoid "Luitpold" in the high dose of 40,000 U. per 100 g cream. By using Hirudoid forte the therapeutically properties unfold their effects by the percutaneous route.

Hirudoid forte cream accelerates the absorption of haematomas and swellings.

Hirudoid forte cream promotes the regression and inhibits the formation of superficial microthrombi.

Hirudoid forte cream has a high affinity to human tissue and is well tolerated.

Composition

Mucopolysaccharide polysulphate

(Organo-heparinoid "Luitpold")

Cream base

0.445 g
ad 100.0 g

Directions for use

Unless otherwise directed by the physician, apply 3-5 cm of cream to the affected area once or twice daily and massage gently until the cream disappears. More cream may be used if necessary. In the case of very painful inflammation, or if there is evidence of thrombosis. Hirudoid forte cream should be applied thickly to the

affected area and covered with gauze or lint. The effect of Hirudoid forte cream can be intensified in these cases by massaging the cream into the skin surrounding the affected area. To soften hard scars, the cream should be massaged vigorously into the scarred area.

Store below 25°C

Presentation

5 g, 7 g, 10 g, 14 g, 20 g, and 40 g of cream

To pierce the aluminium foil at the mouth of the tube, use the spike inside the cap.

Manufactured by:
OLIC (Thailand) Limited
Bangpa-In Industrial Estate, Ayutthaya, Thailand
For: Medinova AG, Switzerland

01883 991009551

 **medinova**
SWITZERLAND